

4th Global Pediatric Ophthalmology Congress

March 07-08, 2019 | Berlin, Germany

Joachim Finger

Roland Consult Stasche & Finger GmbH, Germany

Manfred Stasche

Roland Consult Stasche & Finger GmbH, Germany

Electrophysiology tests for children with RETI-port/scan 21

With device RETI-Port / Scan 21: You will get all programs according to the ISCEV standard. Furthermore it offers the possibility to create own programs, in particular also a short program for pediatric examinations. This is very easy with the stimulators BABYflash E130 and MIniganzfeld I8. With the help of Steady State protocols, one can decide, independently of markers, how the function of the cones and the of the optic nerve are Furthermore there are special fixation pictures for children for the pattern VEP. All data can be exported as.csv file for any scientific statistic analyze.