

Anti-Citrullinated Albumin Antibodies as Biomarker for the Diagnosis of Rheumatoid Arthritis

Vishnupriya R Paturi¹, Ramakrishna R Uppuluri², Lina Gao³, Charles T Roberts¹, Srinivasa R Nagalla^{1*}

¹Diabetomics Incorporated, Hillsboro, Oregon, USA; ²Sri Deepti Rheumatology Center, Hyderabad, India; ³Department of Medicine, Oregon Health and Science University, Portland, Oregon, USA

ABSTRACT

Objective: Rheumatoid Arthritis (RA) is associated with increased levels of citrullinated proteins in synovial fluid and blood. Detection of Anticitrullinated Cyclic Peptide Antibodies (ACPA) is the current standard test for diagnosis of RA in conjunction with clinical symptoms. In addition to ACPA, antibodies to citrullinated human serum albumin (ACA) have been previously reported in individuals with RA. We developed a Point-Of-Care (POC) test for ACA and evaluated its sensitivity and specificity in healthy subjects and subjects with clinically diagnosed RA and other rheumatic conditions and autoimmune disease.

Methods: The ACA POC test is a lateral-flow immunoassay using antihuman IgA/G/M and anti-human serum albumin antibodies for detection of citrullinated serum albumin-ACA complexes. This test was used to determine serum or plasma ACA levels in a South Asian study population comprised of healthy controls (n=484) and subjects with clinically diagnosed RA (n=354) or other rheumatic (n=103) and autoimmune diseases (n=60) and compared to the levels of ACPA and Rheumatoid Factor (RF).

Results: The sensitivity of the ACA POC test for RA was 0.520 and the specificity was 0.994. ACA prevalence in other rheumatoid disease was similar to that of ACPA and less than that of RF. ACA was detected in 12% of RA samples that were negative for ACPA. The combined sensitivity of ACA+ACPA was 0.664 and the combined specificity was 0.845.

Conclusion: The ACA POC test exhibits robust sensitivity and specificity for RA diagnosis in serum or plasma and, in conjunction with ACPA, can enable the rapid and efficient differential diagnosis of RA with increased sensitivity and with comparable specificity.

Keywords: Anticitrullinated cyclic peptide antibodies; Citrullinated human serum albumin antibodies; Diagnosis; Point-of-care; Rheumatoid arthritis; Rheumatoid factor

INTRODUCTION

Rheumatoid Arthritis (RA) is an autoimmune disease with unclear etiology, and is characterized by joint inflammation and underlying bone loss [1]. RA preferentially affects women, and its global prevalence is 1% [1-3]. RA diagnosis is typically accomplished through a combination of clinical assessments that are complemented by determination of the presence of specific biomarkers [4]. RA biomarkers are also potentially important in the identification of at-risk individuals prior to the development of clinically apparent disease in order to institute prophylactic therapies, including new approaches targeting specific inflammatory mediators [5-7]. The principal RA biomarkers are Rheumatoid Factor (RF), which is comprised of antibodies against the Fc domain of IgG, and Anti-Citrullinated Peptide

Antibodies (ACPAs) that recognize proteins containing arginine residues that have been deimided, or citrullinated by members of the Peptidyl Arginine Deimidase (PAD) family [8-11]. This modification is thought to render the citrullinated proteins more immunogenic, resulting in more frequent autoimmune responses [12]. Although RF is detectable in the majority of RA patients, it is also seen in other connective tissue diseases such as Systemic Lupus Erythematosus (SLE) and Sjogren's Syndrome (SS) [13]. Although there is recent evidence that the epitopes recognized by RF and ACPA antibodies may be shared in some cases, ACPAs appear to exhibit more specificity for RA than RF (i.e., up to 80% in serum of RA patients), although ACPAs, like RF, are also detectable in non-RA conditions such as SLE and SS [14-16]. ACPA is the primary biomarker currently employed in RA diagnosis.

A number of studies have characterized a spectrum of citrullinated

Correspondence to: Srinivasa R Nagalla, Diabetomics Incorporated, Hillsboro, Oregon, USA; E-mail: nagallas@diabetomics.com

Received: 04-Jan-2024, Manuscript No. RCR-24-29024; **Editor assigned:** 08-Jan-2024, PreQC No. RCR-24-29024 (PQ); **Reviewed:** 22-Jan-2024, QC No. RCR-24-29024; **Revised:** 29-Jan-2024, Manuscript No. RCR-24-29024 (R); **Published:** 05-Feb-2024, DOI: 10.35841/2161-1149.24.14.378

Citation: Paturi VR, Uppuluri RR, Gao L, Roberts CT, Nagalla SR (2024) Anti-Citrullinated Albumin Antibodies as Biomarker for the Diagnosis of Rheumatoid Arthritis. *Rheumatology (Sunnyvale)*. 14: 378.

Copyright: © 2024 Paturi VR, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

proteins (“citrullinome”) in synovial fluid [17-20] and serum [20] of RA patients. A subset of antibodies against these proteins is what is detected by the standard ACPA test. Citrullinated Human Serum Albumin [HSA] in particular was identified in RA synovial fluid and serum, and Hefton et al subsequently reported the presence of Anti-Citrullinated HSA Antibodies (ACA) in RA serum [18-21]. In that study of 79 RA patients, 38% had ACA levels >2 SD's above healthy controls, and ACA positivity trended higher in females and in patients with joint erosions. In the present study, we extend these observations using a newly developed Point-Of-Care (POC) immunoassay for ACA in serum or plasma and compared its performance to standard ACPA or RF assays in a large East Asian study population.

MATERIALS AND METHODS

During the development of a prototype lateral-flow immunoassay for SARS-CoV-2 antibodies that incorporated HSA detection in the test line, we detected a positive reaction with RA samples. Subsequent investigation determined that the test line positivity reflected the presence of citrullinated HSA-ACA complexes in the RA samples, similar to the earlier finding of Hefton et al [21]. We subsequently designed a lateral-flow immunochromatographic test (described below) to directly assess the levels of ACA in serum or plasma.

Anti-citrullinated HSA antibody lateral-flow POC test

Serum or plasma samples were analyzed using an ACA POC test system (Diabetomics, Inc., Hillsboro, OR). The ACA test system consists of an immunochromatography test strip in a cassette and a portable reader for quantification. Test strips are configured with antibodies against human IgG, IgM, and IgA and human serum albumin. Colloidal gold was used as detection reagent. Streptavidin coupled to colloidal gold and biotin on the control line. 150 µL of 1:80 diluted serum or plasma is added to the test strip and inserted into a proprietary reader that utilizes image analyses for quantification and a QRcode recognition module for test and lot-specific information for calibration and quantification. The ACA concentration is displayed at the end of 15 minutes. Results are reported as arbitrary Units/mL, with a value >20 considered positive (>2 SDs above the range seen in healthy controls). The dynamic range of the ACA assay is 15-4,000 U/mL, and intra- and inter-assay coefficients of variation are 6.2 and 9.8%, respectively.

ACA in subjects with rheumatic and other autoimmune diseases

In a prospective study, serum or plasma samples were collected from 354 subjects with clinically diagnosed RA, 103 with other rheumatic disease (non-RA; 25 SLE, 22 spondyloarthritis, 11 gout, 8 psoriatic arthritis, 7 osteoarthritis, 2 ankylosing spondylitis, 1 polyarthritis, and 27 non-RA rheumatic disease samples of undefined category), and 60 with other autoimmune disease (20 Graves' disease, 10 myasthenia gravis, 20 GADA+ type-1 diabetes, and 10 chronic kidney disease). Samples were tested with the ACA POC device described above as well as with commercial assays for RF and ACPA (RF and anti-CCP IgG assays from Siemens Healthcare Diagnostics; positive values are >12 IU/mL and >5 U/mL, respectively). In a separate cross-sectional study, 484 banked serum (334) or plasma (150) samples from Healthy controls (HC) were tested for ACA. Institutional IRB approval was obtained, the study was registered with regulatory authorities

(CTRI/2020/10/028596), and informed consent was obtained from all study participants.

Statistical analyses

Pair-wise comparisons of ACA levels in the different sample groups employed the Wilcoxon rank sum test followed by Bonferroni adjustment. ACA levels in female vs male RA and non-RA samples were analyzed using the Wilcoxon rank sum test. Age differences in ACA levels were analyzed using one-way ANOVA.

RESULTS

Figure 1A shows ACA levels in HC, RA, and non-RA subjects. Intergroup differences were significant in each pairwise comparison ($p < 0.0001$). At a positive value cut-off of >20 U/mL, 3 of 484 HC subjects (0.6%), 184 of 354 (52%) of RA subjects, and 12 of 103 (12%) non-RA subjects (1 gout, 1 ankylosing spondylitis, 1 osteoarthritis, 1 spondyloarthritis, 6 SLE, and 2 non-RA uncategorized rheumatic disease samples) were positive for ACA. From these data, the calculated sensitivity of ACA for RA was 0.520 and the specificity was 0.996. Figure 1B shows ACA levels in the 184 ACA+ (20 U/mL) RA and 12 ACA+ non-RA subjects. The median ACA levels in the control as shown in Figure 1A and ACA+ RA and ACA+ non-RA as shown in Figure 1B groups were 0.0, 441, and 79 U/mL, respectively. The difference in ACA levels in the ACA+ RA and ACA+ non-RA groups was also significant ($p = 0.019$). In a separate analysis, all 60 samples from patients with Graves' disease, myasthenia gravis, type-1 diabetes, and chronic kidney disease were negative for ACA (data not shown).

We also determined the prevalence of ACPA and RF in the non-RA group. Of the 103 non-RA samples, 8 (8%) were positive for ACPA (SLE, 2 gout, 1 polyarthritis, and 2 uncategorized rheumatic disease) and 24 (23%) were positive for RF (3 gout, 1 polyarthritis, 11 SLE, 4 spondyloarthritis, and 7 uncategorized rheumatic disease). Three of these SLE samples and one of the gout samples were positive for all three biomarkers, but we can't exclude the possibility that these four subjects also had RA. Thus, RF exhibited the most cross-reactivity with non-RA rheumatic disease samples in this sample set. The sensitivity and specificity of ACA, ACPA, and RF in the RA and non-RA sample sets are summarized in Table 1. ACPA and RF values were not available for the HC and other autoimmune disease groups, so ACA specificity in the non-RA group is shown for comparison to the specificity determinations for ACPA and RF using the non-RA group data.

Since ACPA is the primary biomarker currently employed for RA diagnosis, we examined further the relationship between ACA and ACPA status in the RA sample set ($n = 354$). As illustrated in Figure 2, 235 (66%) were positive for either ACPA or ACA. 141 (40%) were positive for both, while 43 (12%) were positive for ACA alone and another 51/354 (14%) were positive for or ACPA alone. Thus, the addition of ACA to ACPA increased the percentage of biomarker-positive samples from 54% to 66% (192 to 235).

We also determined if there was any sex or age difference in ACA levels in the RA and non-RA groups. Figure 3A shows the sex distribution of ACA levels in the RA and non-RA groups. ACA levels were significantly higher in females (Wilcoxon rank sum test p -value=0.0411) in the overall RA group, but not in the non-RA group (Wilcoxon rank sum test p -value=0.205). This is similar to the trend of higher ACA levels in female RA subjects in the US cohort reported previously [21]. Sensitivity of ACA in females

was 0.540 (CI, 0.483, 0.597) and specificity was 0.877 (CI, 0.772, 0.945), similar to the ACA sensitivity and specificity in the total cohort (i.e., 0.520 and 0.883, respectively). ACA sensitivity in males was lower than seen in the overall group (0.372 vs 0.520), while specificity was similar (0.892 vs 0.887). However, as shown in Figure 3B, when the analysis was restricted to ACA+ samples, no

sex difference was apparent (p=0.66 in RA group and p=0.497 in non-RA group).

As shown in Figure 4, there were no significant differences in age between the various biomarker combinations (one-way ANOVA p-value=0.526).





Table 1: Sensitivity and specificity of ACA, ACPA, and RF. Exact binomial confidence limits (95% CI) are shown.

Metric	Estimate	CI
ACA sensitivity in RA	0.52	(0.466, 0.573)
ACA specificity in HC	0.994	(0.982, 0.999)
ACA specificity in non-RA	0.883	(0.805, 0.938)
ACPA sensitivity in RA	0.542	(0.489, 0.595)
ACPA specificity in non-RA	0.922	(0.853, 0.966)
ACPA+ACA sensitivity in RA	0.664	(0.612, 0.713)
ACPA+ACA specificity in non-RA	0.845	(0.76, 0.909)
RF sensitivity in RA	0.689	(0.638, 0.737)
RF specificity in non-RA	0.777	(0.684, 0.853)






Figure 3: ACA levels in female and male subjects in the overall RA and non-RA groups (A): And ACA+RA and non-RA groups (B). **Note:** Sex: (♀) Female; (♂) Male




Figure 4: Age distribution in all biomarker combinations in the RA group.

DISCUSSION

RA is diagnosed and distinguished from other rheumatoid diseases through a combination of clinical assessments and the presence of elevated blood levels of biomarkers, principally ACPA. In the course of development of a rapid POC test for SARS-CoV-2 antibodies, we observed a positive reaction with RA serum samples that was determined to result from the presence of ACA, as previously reported by Hefton et al [21]. The presence of citrullinated HSA in the circulation and the subsequent generation of ACA may be the result of leakage of citrullinated HSA from synovial fluid or direct citrullination of circulating HSA by PAD activity.

PADs are a critical component of Neutrophil Extracellular Trap (NET) formation (NETosis), which is a significant component of RA pathology [22]. Indeed, it was previously shown that the PAD4 isoform implicated in HSA citrullination is released during NETosis from neutrophils in synovial fluid and presumably in the circulation as well, providing a plausible mechanism for generating

citrullinated HSA in the circulation and subsequent generation of ACA [21,23].

These findings prompted us to develop a lateral-flow immunochromatographic POC test to detect ACA. The resulting ACA POC exhibited excellent sensitivity (0.520). And specificity (0.996) in a large East Asian sample set comprised of healthy controls, clinically diagnosed RA cases, and other rheumatoid and autoimmune diseases. ACA levels detected with this test were present in a small percentage (14%) of non-RA rheumatic disease samples, vs 8 and 23% for ACPA and RF, respectively. Importantly, a significant number of RA samples were positive for ACA but not for ACPA. In this sample set, the addition of ACA to ACPA increased the detection of RA from 54 to 66%, thereby increasing overall sensitivity. The addition of an ACPA test line to the current ACA POC would provide a convenient and comprehensive assessment of RA status that would be particularly useful in low-resource environments. Wang et al. have recently described the use of the POC ACA test described here in both serum and plasma

samples from a US cohort, and also found that ACA was highly specific for RA [24].

Our study does have a limitation in that we did not have access to a detailed clinical characterization of disease severity or treatment status and could not, therefore, assess any associations of ACA positivity with these parameters in this patient population. However, Wang et al, using the POC test described herein in a US patient population, recently reported that ACPA+/ACA+ subjects in their study exhibited a trend toward a higher proportion of joint erosions and that ACPA-/ACA+ subjects did not include any TNF inhibitor non-responders *vs* those in the ACPA-/ACA-, ACPA+/ACA-, and ACPA+/ACA+ groups and were more likely to be treated with a Tumor Necrosis Factor (TNF) inhibitor or methotrexate alone.

AUTHOR CONTRIBUTIONS

VPR, CTR, and SRN designed the study, analyzed data, and wrote the manuscript. VPR and RRU generated data. LG analyzed data.

Funding

This study was funded by Diabetomics, Inc.

Conflict of interest

VRP and SRN are employees of, and VRP, CTR, and SRN are shareholders in Diabetomics, Inc., which funded this study.

REFERENCES

1. Scott DL, Wolfe F, Huizinga TW. Rheumatoid arthritis. *Lancet*. 2010;376:1094-1098.
2. Pedersen JK, Kjær NK, Svendsen AJ, Hørslev-Petersen K. Incidence of rheumatoid arthritis from 1995 to 2001: Impact of ascertainment from multiple sources. *Rheumatol Int*. 2009;29(4):411-415.
3. Alamanos Y, Drosos AA. Epidemiology of adult rheumatoid arthritis. *Autoimmun Rev*. 2005;4(3):130-136.
4. Aletaha D, Neogi T, Silman AJ, Funovits J, Felson DT, Bingham III CO, et al. 2010 rheumatoid arthritis classification criteria: An American College of Rheumatology/European League Against Rheumatism collaborative initiative. *Arthritis Rheum*. 2010;62(9):2569-2581.
5. Deane KD, Holers VM. Rheumatoid arthritis pathogenesis, prediction, and prevention: An emerging paradigm shift. *Arthritis Rheumatol*. 2021;73(2):181-193.
6. Haville S, Deane KD. Pre-RA: Can early diagnosis lead to prevention?. *Best Pract Res Clin Rheumatol*. 2022;36(1):101737.
7. Senolt L. Emerging therapies in rheumatoid arthritis: Focus on monoclonal antibodies. *F1000Res*. 2019;8:F1000 Faculty Rev-1549.
8. Rönnelid J, Turesson C, Kastbom A. Autoantibodies in rheumatoid arthritis—laboratory and clinical perspectives. *Front Immunol*. 2021;12:685312.
9. Schellekens GA, De Jong BA, van den Hoogen FH, van de Putte LB, van Venrooij WJ. Citrulline is an essential constituent of antigenic determinants recognized by rheumatoid arthritis-specific autoantibodies. *J Clin Invest*. 1998;101(1):273-281.
10. Vossenaar ER, Zendman AJ, van Venrooij WJ, Pruijn GJ. PAD, a growing family of citrullinating enzymes: genes, features and involvement in disease. *Bioessays*. 2003;25(11):1106-1118.
11. Tilwala R, Thompson PR. Peptidyl arginine deiminases: Detection and functional analysis of protein citrullination. *Curr Opin Struct Biol*. 2019;59:205-215.
12. Nguyen H, James EA. Immune recognition of citrullinated epitopes. *Immunology*. 2016;149(2):131-138.
13. Nell V, Machold KP, Stamm TA, Eberl G, Heinzl H, Uffmann M, et al. Autoantibody profiling as early diagnostic and prognostic tool for rheumatoid arthritis. *Ann Rheum Dis*. 2005;64(12):1731-1736.
14. Trela M, Perera S, Sheeran T, Rylance P, Nelson PN, Attridge K, et al. Citrullination facilitates cross-reactivity of rheumatoid factor with non-IgG1 Fc epitopes in rheumatoid arthritis. *Sci Rep*. 2019;9(1):12068.
15. Trier NH, Holm BE, Hansen PR, Slot O, Loch H, Houen G, et al. Specificity of anti-citrullinated protein antibodies in rheumatoid arthritis. *Antibodies (Basel)*. 2019;8(2):37.
16. Payet J, Gouvestre C, Bialé L, Avouac J, Wipff J, Job-Deslandre C, et al. Anticyclic citrullinated peptide antibodies in rheumatoid and nonrheumatoid rheumatic disorders: Experience with 1162 patients. *J Rheumatol*. 2014;41(12):2395-2402.
17. van Beers JJ, Schwarte CM, Stammen-Vogelzangs J, Oosterink E, Božič B, Pruijn GJ, et al. The rheumatoid arthritis synovial fluid citrullinome reveals novel citrullinated epitopes in apolipoprotein E, myeloid nuclear differentiation antigen, and β -actin. *Arthritis Rheum*. 2013;65(1):69-80.
18. Bennike T, Lauridsen K, Olesen MK, Andersen V, Birkelund S, Stensballe A, et al. Optimizing the identification of citrullinated peptides by mass spectrometry: Utilizing the inability of trypsin to cleave after citrullinated amino acids. *J Proteomics Bioinform*. 2013;6:288-295.
19. Tuttüren AE, Fleckenstein B, de Souza GA. Assessing the citrullinome in rheumatoid arthritis synovial fluid with and without enrichment of citrullinated peptides. *J Proteome Res*. 2014;13(6):2867-2873.
20. Tilwala R, Nguyen SH, Maurais AJ, Nemmara VV, Nagar M, Salinger AJ, et al. The rheumatoid arthritis-associated citrullinome. *Cell Chem Biol*. 2018;25(6):691-704.
21. Hefton A, Liang SY, Ni K, Carter V, Ukadike K, Lood C, et al. Autoantibodies against citrullinated serum albumin in patients with rheumatoid arthritis. *J Transl Autoimmun*. 2019;2:100023.
22. Song W, Ye J, Pan N, Tan C, Herrmann M. Neutrophil extracellular traps tied to rheumatoid arthritis: Points to ponder. *Front Immunol*. 2021;11:578129.
23. Spengler J, Lugonja B, Jimmy Ytterberg A, Zubarev RA, Creese AJ, Pearson MJ, et al. Release of active peptidyl arginine deiminases by neutrophils can explain production of extracellular citrullinated autoantigens in rheumatoid arthritis synovial fluid. *Arthritis Rheumatol*. 2015;67(12):3135-3145.
24. Wang X, Laine A, Chung S, Mustelin T. Diagnostic and prognostic utility of anti-modified albumin autoantibodies in rheumatoid arthritis. *Arch Rheum Arth Res*. 2022;2(2):ARAR.MS.ID.000535.