

A Quick and Simple HPLC Method for the Quantification of the Immunosuppressor 6-Mercaptopurine in Human Plasma

Claudia Schildböck*, Stephan Harm and Jens Hartmann

Department for Biomedical Research, Danube University, Krems, Austria

Abstract

Thiopurine drugs, including Azathioprine and 6-Mercaptopurine (6-MP), are used in clinical practice for more than 50 years. Their bioavailability is strongly patient-dependent. Therefore, closely monitoring of the plasma levels during therapy is recommended. Here, we present a quick and simple High Performance Liquid Chromatography (HPLC) based method for the quantification of 6-MP in human plasma. Samples were deproteinized with methanol to obtain protein concentrations below 2 mg/dL. The method offers a recovery higher than 90% and is precise and linear with a high correlation coefficient and a level of detection <100 ng/mL. The runtime of five minutes offers short analysis latency. This HPLC based method can be used for monitoring the 6-MP serum level in human blood.

Keywords: 6-Mercaptopurine; HPLC; Immunosuppressor

Introduction

6-Mercaptopurine was developed in the 1950s as an anticancer agent [1]. Some years later, it could be shown that 6-MP has also an immunosuppressive effect [2]. From the early 1950s to the early 1980s, Azathioprine, the prodrug of 6-MP, was one of the main immunosuppressive drugs in renal transplantation. Due to low costs, Azathioprine and 6-MP are still used in the treatment of acute lymphoblastic leukemia and inflammatory bowel diseases such as Crohn's disease or Ulcerative Colitis.

Orally applied Azathioprine is non-enzymatically converted to 6-MP, which is able to pass the cell membrane. There are three intracellular pathways which compete with each other depending on the availability of various activating and inactivating enzymes or cofactors. One of these pathways converts 6-MP to metabolic active 6-thioguanine nucleotides which are incorporated into DNA and RNA instead of adenine and guanine leading to suppressed cell proliferation, mainly of lymphocytes. The other two pathways end up building inactive metabolites [3,4] (Figure 1).

Since the amounts of enzymes and cofactors which play a role in the metabolism of 6-MP show a wide interindividual variation, the response to Azathioprine treatment is unpredictable [4,5]. In order to maintain a constant serum level of 6-MP, analytical monitoring of the drug could be of great advantage. Currently, there are no test kits or recommended reference methods for the clinical monitoring of 6-MP commercially available. Therefore, numerous HPLC methods for the determination of 6-MP in human erythrocytes, plasma and pharmaceutical formulations have been published. Most of these procedures combine the detection of 6-MP with its metabolites. Other methods show complicated and time-consuming sample preparation procedures, such as HPLC-tandem mass spectroscopy and extraction with chloroform combined with the use of special and expensive equipment and software [6-8].

The objectives of this study were to create a quick and simple HPLC method with good recovery, selectivity and accuracy for the quantification of 6-MP in human plasma well as an optimized sample preparation which minimizes the protein content in the sample in order to gain a long column life-time.

Materials and Methods

Instrumentation

The HPLC device was a Shimadzu prominence system (Kyoto, Japan) consisting of a DGU-20A-5R degassing unit, a LC-20AD pump unit, a SIL-20AC-HT autosampler, a CTO-20AC column oven and an SPD-20AV UV/VIS detector. For data processing, LabSolutions software version 5.71 was used. Short thread clear 1.5 mL glass HPLC vials with ultra clean seal 9 mm PP caps with center hole were purchased from Fisher Scientific (Hampton, New Hampshire, USA). The analytical column was a Nukleosil 100-5C18, 4.6 × 150 mm purchased from Macherey-Nagel (Duren, Germany) which was combined with a 4.6 × 5 mm Onyx Monolithic C18 guard column (Phenomenex, Torrance, California, USA) to reduce contaminations.

Chemicals and reagents

6-MP and 5-sulfosalicylic acid were obtained from Sigma Aldrich (St. Louis, Missouri, USA), HPLC grade methanol was purchased from Carl Roth (Karlsruhe, Germany) and distilled water was obtained from Fresenius Kabi (Graz, Austria).

For each analytical day a 6-MP 200 µg/mL stock solution was prepared in methanol.

Protein binding of 6-mercaptopurine

Since only the amount of freely available and unbound fraction is chromatographically detectable, the protein binding of 6-mercaptopurine was determined. For this purpose, plasma was spiked with 6-mercaptopurine, incubated for 30 minutes and centrifuged for 20 minutes at 15,000 g through a centrifugal concentrator with a molecular cut off of 30 kDa (Spin-X UF 500 µL, Corning, New York, USA). The levels of 6-MP were determined in spiked plasma and in the

*Corresponding author: Claudia Schildböck, Department for Biomedical Research, Danube University, Krems, Austria, Tel: +4327328932619; E-mail: claudia.schildboeck@donau-uni.ac.at

Received April 30, 2019; Accepted May 16, 2019; Published May 23, 2019

Citation: Schildböck C, Harm S, Hartmann J (2019) A Quick and Simple HPLC Method for the Quantification of the Immunosuppressor 6-Mercaptopurine in Human Plasma. J Chromatogr Sep Tech 10: 420.

Copyright: ©2019 Schildböck C, et al. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

filtrate. The percentage of protein bound fraction was calculated using the formula:

$$\text{Protein binding [\%]} = \frac{C_{\text{filtrate}}}{C_{\text{plasma}}} \times 100 \quad (1)$$

C_{plasma} is the determined concentration of 6-MP and C_{filtrate} is the measured drug level in filtrate.

6-MP extraction from plasma samples

In order to dissolve the protein-bound part of 6-MP on the one hand and not to load the column with proteins on the other hand to prolong the column's live time, protein precipitation was performed for all plasma samples. Although many possibilities for the extraction of drugs from human plasma have been described in the literature [9], only methanol and 5-sulfosalicylic acid, two commonly used reagents for protein precipitation [10], were compared regarding their deproteinization efficiency. Methanol can be mixed with water in any ratio and its low dielectric constant reduces the solvation power of the aqueous solution for dissolved proteins. At temperatures above 10°C a considerable denaturation of the proteins can occur; therefore,

precipitation should be carried out at a temperature below 0°C [11]. To extract 6-MP from human plasma, four different approaches were conducted (Table 1).

Two different solutions with 0.5 g/mL 5-sulfosalicylic acid solutions in water and 50% methanol were prepared. Precipitation was carried out by mixing 100 µL of the prepared solution with 1 mL plasma, incubated at 4°C for 20 minutes and centrifuged at 15,000 g.

With methanol two different extraction procedures were tested. According to the recommended extraction method by Naik et al. [12], 100 µL plasma were mixed with 1 mL methanol to water (1:1), vortexed for 5 minutes and centrifuged at 3000 g for 15 minutes at room temperature. In comparison, precipitation at low temperature was tested where 100 µL plasma were mixed with 500 µL iced methanol (-20°C), vortexed and incubated for 20 minutes at -80°C followed by centrifugation at 15,000 g for 5 minutes (Table 1). Protein quantification of the supernatant of all precipitated plasma samples was carried out with a Hitachi cobas c311 automated analyzer with the corresponding test kit, both purchased from Roche (Basel, Switzerland). For the HPLC method described below, protein precipitation was performed with iced




Figure 1: Metabolic pathway for Azathioprine: The initial step in 6-mercaptopurine (6-MP) nucleotide metabolite formation is catalysed by hypoxanthine phosphoribosyltransferase (HPRT). Further metabolism of mercaptopurine nucleotide results in the eventual formation of mercaptopurine nucleotide. Oxidation catalysed by xanthine oxidase, and methylation catalysed by thiopurine methyltransferase, compete with thioinosine triphosphate and thioidoxyguanosine triphosphate formation. Thioidoxyguanosine metabolites of 6-MP act as purine antagonists and inhibit DNA, RNA, and protein synthesis, inducing cytotoxicity and immunosuppression.

methanol in the ratio plasma to methanol 1:5, as this sample preparation yielded the lowest protein concentration in the supernatant (Table 2).

Chromatographic conditions

Methanol and water were used as mobile phase. The water content in the solvent is necessary for the retention of organic substances on a reversed phase column. The methanol content controls the retention time and ensures good transparency in the UV/VIS detector [12]. Using an isocratic mixing ratio of the mobile phase, it was not possible to separate the 6-MP peak from other plasma compounds. Therefore, the mobile phase gradient was empirically optimized to get a clear separation. The column temperature was maintained at 25°C during

chromatography. The sample injection volume was 20 µL, and the run time per sample was 5 minutes using a flow rate of 0.8 mL/min. The composition of the mobile phase started with a methanol to water ratio of 9:1 and changed following a linear gradient with the ratio 7:3 after 1.7 minutes. The eluted 6-MP peaks were detected via UV/VIS detector at 325 nm [13]. Peak area was used to determine 6-MP concentrations in plasma samples.

Preparation of calibration curves



A 200 µg/mL 6-MP standard solution was prepared and diluted in mobile phase (water to methanol 1:9) to obtain 6-MP concentrations of 1000, 750, 500, 250, 100 and 0 ng/mL. Before the standards were injected into the HPLC system, they were diluted with a ratio (plasma to methanol) of 1:5 to compensate the dilution of the plasma samples during the protein precipitation step. For the preparation of the calibration curve in human plasma, plasma samples were spiked with 6-MP stock solution to obtain final concentrations of 1000, 750, 500, 250, 100 and 0 ng/mL. For this procedure, the 6-MP stock solution (200 µg/mL in methanol) was diluted with physiological sodium chloride solution to get concentrations of 10, 7.5, 5.0, 2.5, 1 and 0 µg/mL. Each solution was then diluted 1:10 with plasma to obtain the desired standard concentrations. This procedure ensured that each standard contained the same concentration of plasma, which enabled a more precise and reproducible peak integration. Protein precipitation was then performed as described above and samples were injected into the HPLC system after centrifugation. For both calibration curves, linear plots of peak area against 6-MP concentration were created.

Statistical analysis

All statistical calculations were carried out in Microsoft Excel 2016 for Windows.

Validation parameters

Recovery: Recovery was determined by comparison of the 6-MP peak areas spiked in mobile phase with human plasma. Additionally, four quality control (QC) samples were prepared at the concentrations of 800, 400, 200 and 100 ng/mL 6-MP in human plasma. All samples



Precipitation solution	Ratio [v/v] precipitation solution: Plasma	Incubation temperature
5-Sulfosalicylic acid in water	1:10	4°C
5-Sulfosalicylic acid in water: methanol 1:1	1:10	4°C
Methanol: water 1:1	10:1	20°C
Iced methanol	5:1	-80°C

Table 1: Four different protein precipitation protocols were compared regarding protein content in the supernatant.

Precipitation solution	Total protein concentration [mg/dL]
5-Sulfosalicylic acid in water	101.5
5-Sulfosalicylic acid in water: methanol 1:1	83.8
Methanol: water 1:1	390.9
Iced methanol	1.6
Precipitation solution	Total protein concentration [mg/dL]
5-Sulfosalicylic acid in water	101.5
5-Sulfosalicylic acid in water: methanol 1:1	83.8
Methanol: water 1:1	390.9
Iced methanol	1.6

Table 2: Comparison of the protein concentration in the supernatant for different protein precipitation methods for 6-MP spiked plasma.

Nominal conc. [ng/mL]	Interday (n=4)				Intraday (n=5)			
	Precision		Recovery	Accuracy	Precision		Recovery	Accuracy
	Mean [ng/mL]	RSD [%]	[%]	MRE [%]	Mean [ng/mL]	RSD [%]	[%]	MRE [%]
1000 (STD)	900.9	11.9	90.1	5.5	900.9	8.2	90.1	7.9
750 (STD)	748.3	13.9	99.8	2.2	748.3	7.6	99.8	1.3
500 (STD)	509.1	8.4	101.6	-3.6	512.0	5.7	102.4	-3.5
250 (STD)	285.2	7.1	114.4	-11.6	293.9	4.7	117.6	-13.0
100 (STD)	98.1	22.2	98.1	-5.5	88.9	19.0	88.9	-7.8
800 (QC)	725.4	15.5	90.7	1.3	725.4	5.3	90.7	4.2
400 (QC)	372.3	6.2	93.1	-7.2	372.3	4.0	93.1	-7.2
200 (QC)	245.9	11.1	123.0	-13.8	245.9	3.9	123.0	-13.8
100 (QC)	104.8	20.7	104.8	17.5	90.1	8.6	90.7	10.3

RSD [%]: Relative Standard Deviation (=Coefficient of variation), MRE [%]: Mean Relative Error, STD: Standard Samples, QC: Quality Control Samples.

Table 3: Summary of the statistic parameters for the analysis of 6-MP.

Variables	RSD%			
	Plasma samples	4°C	- 20°C	- 80°C
Freeze thaw cycles (n=3)	QC 100 ng/mL	-	26	50
	QC 200 ng/mL	-	5	4
	QC 400 ng/mL	-	2	2
	QC 800 ng/mL	-	3	5
Variables	Recovery %			
	Plasma samples	4°C	- 20°C	- 80°C
28-day stability (n=3)	QC 100 ng/mL	3972	296	176
	QC 200 ng/mL	818	181	153
	QC 400 ng/mL	174	129	67
	QC 800 ng/mL	112	107	108

Table 4: RSD values in % of QC samples after 3 freeze thaw cycles and recovery in % of QC samples after 28-day storage at different temperatures.

were measured five times daily (intraday) on four consecutive days (interday).

Accuracy: Accuracy was determined by comparing the measured concentrations in plasma with the spiked concentration. It was expressed as mean relative error in percent which was calculated with the formula:

$$rE [\%] = \frac{X_m - X_s}{X_s} \times 100 \quad (2)$$

rE is the relative error, X_m is the measured concentration and X_s is the spiked concentration. All samples were measured five times a day on four consecutive days.

Precision: The precision was evaluated by performing replicate analyses of QC and standard samples in human plasma. The precision was calculated as the relative standard deviation (RSD%) of interday and intraday measurements. For intraday precision, all measurements were repeated five times on the same day. For interday precision, all measurements were replicated on 4 consecutive days. Then RSD% values were calculated for intraday and interday measurements.

Stability: The stability of the plasma samples containing 800, 400, 200, 100 ng/mL 6-MP were assessed in triplicates at 4°C, -20°C, -80°C for 28 days.

Freeze-thaw stability of these samples was evaluated in triplicates after three complete freeze-thaw cycles from -20°C to room temperature and from -80°C to room temperature on successive days.

Level of detection and level of quantification: According to the ICH guidelines, the Level of Detection (LOD) and the Level of Quantification (LOQ) were assessed based on the Standard Deviation (SD) of the nominal concentrations of the calibration curve, whereas $LOD=3 \times SD$ and $LOQ=10 \times SD$ [13].

Results and Discussion

Protein binding of 6-mercaptopurine

Binding studies of 6-MP to human serum albumin show that between 20% and 90% of the drug is irreversibly bound to albumin [14]. In our current study the determined value of bound 6-MP to plasma proteins was $33.8 \pm 6.3\%$. This fact shows the importance of a good extraction procedure of 6-MP in human plasma (Table 2).

6-MP extraction from plasma samples

Plasma samples which were precipitated 1:10 with methanol to water (1:1) and then centrifuged for 15 minutes at room temperature showed the highest protein concentration left in the supernatant. Both approaches with 5-sulfosalicylic acid showed a high protein concentration left in the supernatant as well. Furthermore, sulfosalicylic acid caused an additional peak which overlaid the 6-MP peak. Plasma which was precipitated with iced methanol had the lowest protein concentration in the supernatant and did not show an interfering peak in the chromatogram (Table 2 and Figure 2). This extraction method is also recommended in the literature for the sample preparation of human plasma samples for other drugs [15].

Plasma proteins were most efficiently removed by methanol precipitation at a ratio of 1:5 v/v for plasma to methanol. Furthermore, since methanol is also part of the mobile phase, there are no disturbing peaks in comparison to sulfosalicylic acid. Therefore, methanol precipitation not only prolongs the lifetime of the HPLC column, but also simplifies data assessment. The following procedure was chosen for the development of this HPLC method: 100 μ L plasma samples were precipitated with 500 μ L iced methanol, vortexed and frozen at -80°C for 20 minutes. After centrifugation at 15,000 g for 5 minutes, 20 μ L of the supernatant were injected. The 6-MP peak was detected at a retention time of 2.25 ± 0.3 minutes (Figure 2).

Calibration curves (linearity)

The calibration curves for quantification of 6-MP in mobile phase as well as in human plasma were developed based on six different concentrations (0, 100, 250, 500, 750 and 1000 ng/mL) using a linear trend line with the formula:

$$y = k \times x + d$$

which is depicted as peak area (x-axis) versus 6-MP concentration (y-axis) in ng/mL. Both curves show an excellent linearity with a least square fit (R^2) of > 0.99 (Figure 3).

Statistical parameters

For the validation of the method, the statistic parameters recovery, accuracy and precision were determined within an analysis day (intraday) and day by day (interday).

Recovery

The concentrations of 6-MP were calculated with the standard curve for mobile phase. The interday recovery was 90%, 100%, 102%, 114% and 98% with a concentration of 1000 ng/mL, 750 ng/mL, 500 ng/mL, 250 ng/mL and 100 ng/mL, respectively (Table 3). Deviation was less than 20% for all concentrations [16]. These findings confirm the efficient protein precipitation for the plasma samples.

Accuracy

In order to compare the accuracy of the measurements, the mean relative error (MRE) was determined. The average MRE for all nominal 6-MP concentrations shown in Table 3 was -2.0%, while the largest MRE was 17.5%.

Precision

Standard plasma samples showed a mean interday concentration of 901 ± 107 ng/mL, 748 ± 127 ng/mL, 509 ± 42 ng/mL, 285 ± 25 ng/mL and 98 ± 42 ng/mL and a mean intraday concentration of 901 ± 74 ng/mL, 748 ± 57 ng/mL, 512 ± 29 ng/mL, 294 ± 19 ng/mL and 89 ± 22 ng/mL. Mean standard deviation was less than 20% in all plasma samples except samples containing 100 ng/mL when measured on different days (Table 3).

Quality control samples had a mean interday concentration of 725 ± 136 ng/mL, 372 ± 132 ng/mL, 246 ± 26 ng/mL and 105 ± 96 ng/mL and a mean intraday concentration of 725 ± 50 ng/mL, 372 ± 17 ng/mL, 246 ± 9 ng/mL and 90 ± 28 ng/mL. Here also the samples containing 100 ng/mL measured on different days showed a higher mean standard deviation than 20%, all others were within the acceptable 20% [16]. All concentrations were within the range of acceptable deviation, only concentrations of 100 ng/mL could not be quantified satisfyingly (RSD% >20%).

Stability

Freeze thaw cycles showed a relative standard deviation of less than 5% with no difference if the samples were frozen at -20 or -80°C for plasma samples ranging from 200 to 800 ng/mL 6-MP. Higher RSD% was only observed in plasma samples with 100 ng/mL (Table 4). The stability test showed that plasma samples with low 6-MP concentrations are not able to be stored for 28 days. In addition, it could be shown that recovery at a storage temperature of -80°C is better than at -20°C. The plasma samples should not be stored in the refrigerator at 4°C (Table 4). The results of the stability tests indicate that plasma samples should be extracted and analyzed immediately after sampling.

Level of detection and level of quantification

Although the detection limit (LOD) of the method presented here is below 100 ng/mL, we consider the method to be valid for the quantification (LOQ) of 6-MP for concentrations of at least 100-200 ng/mL. For the quantification of lower 6-MP serum levels, more sensitive methods based on mass spectroscopy [17] should be considered.

Conclusion

Although the reversed phase HPLC based method developed within this study is not as sensitive as LC-MS based methods, it is adequate for the determination of 6-MP serum levels in human blood, especially since the sample preparation is cheap and easy to handle, and the short runtime of 5 minutes leads to quick results. However, for a better estimation of the *in vivo* 6-MP therapy success, the active metabolite 6-thioguanine in erythrocytes can be performed additionally.

References

1. Hitchings GH, Elion GB (1954) The chemistry and biochemistry of purine analogs. Ann NY Acad Sci 60: 195-199.
2. Schwartz R, Dameshek W (1960) The effects of 6-mercaptopurine on homograft reactions. J Clin Invest 39: 952-958.
3. Reynolds NJ, Sinha A, Elias MS, Meggitt SJ (2016) Translating translation into patient benefit for atopic eczema. Br J Dermatol 175 Suppl 2: 8-12.
4. Supandi S, Harahap Y, Harmita H, Andalusia R (2018) Quantification of 6-mercaptopurine and its metabolites in patients with acute lymphatic leukemia using dried blood spots and UPLC-MS/MS. Sci Pharm 86: 1-18.
5. Schmiegelow K, Kriegbaum NJ (1993) 6-Thioguanine nucleotide accumulation in erythrocytes during azathioprine treatment for systemic connective tissue diseases: A possible index for monitoring treatment. Ann Rheum Dis 52: 152-154.
6. Mohammad-Hossein S, Mohammad-Yaser K, Kaveh A, Abdolhossein N, Davoud A, et al. (2011) Simultaneous determination of 6-mercaptopurine and its oxidative metabolites in synthetic solutions and human plasma using spectrophotometric multivariate calibration methods. Bioimpacts 1: 53-62.
7. Hawwa AF, Millership JS, Collier PS, McElnay JC (2009) Development and validation of an HPLC method for the rapid and simultaneous determination of 6-mercaptopurine and four of its metabolites in plasma and red blood cells. J Pharm Biomed Anal 49: 401-409.
8. Kirchherr H, Kuhn-Velten WN (2005) HPLC-tandem mass spectrometric method for rapid quantification of dimethylarginines in human plasma. Clin Chem 51: 249-52.
9. Soda M, Shibata Y, Yasue M, Fujimura M, Takahashi H, et al. (2015) Simple HPLC Method for the Determination of Caspofungin in Human Plasma. Clin Pharmacol Biopharm 4.
10. Block WD, Markovs ME, Steele BF (1966) Comparison between free amino acid levels in plasma deproteinated with picric acid and with sulfosalicylic acid. Proc Soc Exp Biol Med 122: 1089-1091.
11. Alshammari TM, Ali Al-Hassan A, Hadda TB, Mohamad A (2015) Comparison of different serum sample extraction methods and their suitability for mass spectrometry analysis. Saudi Pharm J 23: 689-697.
12. Naik KM, Sharanappa T, Nandibewoor P (2013) RP-HPLC method for the estimation of 6-mercaptopurine in spiked human plasma and pharmaceutical formulations. Journal of Analytical Chemistry 68: 1085-1088.
13. Baber N (1994) International Conference on harmonisation of technical requirements for registration of pharmaceuticals for human use (ICH). Br J Clin Pharmacol 37: 401-404.
14. Sjöholm I, Stjerna B (1981) Binding of drugs to human serum albumin XVII: Irreversible binding of mercaptopurine to human serum proteins. J Pharm Sci 70: 1290-1291.
15. Zhang X, Zhai S, Duan J, Yang L (2016) A validated HPLC-MS/MS method for the quantification of caderofloxacin in human plasma and its application to a clinical pharmacokinetic study. Biomed Chromatogr 30: 126-130.

16. UNODC (2009) Guidance for the validation of analytical methodology and calibration of equipment used for testing of illicit drugs in seized materials and biological specimens. New York: United Nations Office on Drugs and Crime. p: 76.
17. Cuffari CS, Hunt S, Bayless TM (2000) Enhanced bioavailability of azathioprine compared to 6-mercaptopurine therapy in inflammatory bowel disease: Correlation with treatment efficacy. *Aliment Pharmacol Ther* 14: 1009-1014.