

3rd International Conference and Exhibition on **Cosmetology & Trichology**

July 21-23, 2014 Hampton Inn Tropicana, Las Vegas, USA

Fill/Traction to Fill/Stim

Pier Antonio Bacci
Roma International Estetica, Italy

The process of aging of the face is caused by a series of mechanical, physics, and physiological chemical alterations that reduce the normal metabolic activities and the vascular functions. So we have those alterations and chronic degenerations characterized by cutaneous excess, reduction and alteration off at tissue, prolapse of the skin. According to the type and of the degree of ageing of the face different treatments can be used, such as surgical or dermatologic cosmetic treatments, on respect of different aesthetic pathologies. In advanced cases, we have traditional lifting but today we can use different mini invasive strategies too. In patients without a true cutaneous excess and tissue prolapse and where the aesthetic pathology is characterized by a reduction of the "structure" and of the "tropism" of the skin, traditional lifting is not always suitable and is not the first indicated treatment, even if it can be sometimes of help. These cases, as a rule, represent the large majority of the patients of 35 to 50 years old. The name "Fill-Traction" means an integrate protocol of treatment for a new youth of the face using in the same session fillers, botulin toxin and absorbable unidirectional barbed threads (Medical Fill-Traction). After 10 years of experience this integrate methodology provides us good results improving the quality of the skin, reducing wrinkles and folds, and harmonizing external aspect. Today we can use a not invasive methodology by PDO micro threads (normal, cog or screws) providing us an interesting stimulation of the tissues, thanks to the improvement of the vascularization and cellular activity. This strategy of bio stimulation may be improved by a periodical use of external oxygen propulsion in association to an external not ablative laser power by a diode laser at 1470 nanometers. This way provides us the stimulation of the mitochondrial activity and the cellular metabolism, the use of the PDO micro threads, inserted as a "quilting way", provides us an interesting and topical rejuvenation of the tissues

baccipa@ntc.it