

Final Implanting and Test for Errection Penile Implant Prosthesis

Sabri Dib*

Chirurgien Urologue, Cité 'Le Sahel' - Bt A18, Kolea, W. Tipaza, Algeria

*Corresponding author: Sabri Dib, Chirurgien Urologue, Cite 'Le Sahel'-Bt A18, Kolea, W. Tipaza, Algeria, Tel: +213 20 530 377; Fax: +213 24 483 085; E-mail: dib96dz@yahoo.fr

Received date: Mar 10, 2015, Accepted date: Mar 12, 2015, Published date: Mar 16, 2015

Copyright: © 2015 Dib S. This is an open-access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Medical Image



Figure 1: Final implanting and test for erection penile implant prosthesis.

In the end; all parts of the device are implanted inside the body after sutured the corpora using absorbable suture.

The patient simply pump it up for an erection and press the release valve above the pump to deflate.

We have to signal the importance of the pressure dressing in the first 24 postoperative hours.

The patient usually leaves the hospital in 24 hours.

When the patient returned home and edema secondary to intervention disappeared ,we can ask the patient to pull down the pump located in the scrotum to adjust it.

After a period of 3 to 6 weeks; the patient can start operating the unit for the first time (inflating and deflating the prosthesis).

Four to six weeks after surgery the patient can begin to use the prosthesis for Sexual Relations.